




ONE-PIECE SERIES | MBI NCT™

RECOMMENDED STRAIGHT DRILL PROTOCOL

Drill Diameter [mm]		Ø1.2	Ø1.5	Ø2.0
Drill Speed [RPM]		1200-1500	1200-1500	900-1200
IMPLANT DIAMETER	Ø2.0	Soft Bone	2/3	↓
		Hard Bone	↓	
	Ø2.4	Soft Bone	↓	→ 2/3 ↓
		Hard Bone		↓
	Ø2.9	Soft Bone		↓
		Hard Bone	↓	→ 2/3 ↓


 Drill to mark osteotomy site

 Drill osteotomy to implant

 Drill osteotomy partially according to implant

The recommended drill protocol procedure should not replace the dentist's/surgeon's judgment. The implants may be loaded for immediate function when good primary stability (above 35 Ncm) has been achieved and with appropriate occlusal loading.

MBI & MBI NC COMPONENTS



The diagram displays the following components:

- NM-T3201**: Metal Housing (Dimensions: $\phi 4.3$, 2.8)
- NM-T3202**: Extra Soft Nylon Cap
- NM-T3203**: Soft Nylon Cap (Dimensions: $\phi 3.4$, 2.3, MICRO)
- NMCD2312**: Drill $\phi 1.2\text{mm}$
- NMCD2315**: Drill $\phi 1.5\text{mm}$
- NM-X1520**: Driver $\phi 1.7\text{mm}$
- NM-X2415**: Driver $\phi 1.7\text{mm}$
- NM-T6124**: Analog
- NM-T4410**: Transfer