


CHALLENGE SERIES | PTERYCORE™



BONE TYPES	Pterygoid Region
PROSTHETICS PLATFORM	Internal hex
DESIGN FEATURES	<ul style="list-style-type: none"> • Large cutting surface area • Tapered thread and tapered core body • Wide and sharp threads • Machined surface and reduced diameter "neck"
CLINICAL BENEFITS	<ul style="list-style-type: none"> • For pterygomaxillary region (posterior region of the atrophic maxilla) • Self tapping • High primary stability • Reduces adherence of perio-pathogens thus minimizes the chances for inflammation • Reduced pressure on "neck" portion



ORDERING INFORMATION

Ø D (mm)	Ø D1 (mm)	L (mm)	Ref. No
4.2	2.5	18	NM-M4218
		20	NM-M4220
		22.5	NM-M4222
		25	NM-M4225

Cover Screw Included with all implants



NM-S5023

RECOMMENDED STRAIGHT DRILL PROTOCOL

IMPLANT DIAMETER	Drill Diameter [mm]	Osteotome	Ø2.0	Ø2.8
	Drill Speed [RPM]		1200-900	1000-800
Ø4.2	All Bone Types	1/2 → → →		
↓	Drill osteotomy to implant			
	Drill osteotomy partially according to implant			

Prior to the use of Pteryfit/Pterycore implants in the Pterygoid regions, additional training is recommended.