


CEMENT RETAINED & SCREW RESTORATIONS | ESTHETIC SCREW ABUTMENTS





DESIGN FEATURES	<ul style="list-style-type: none"> • A range of cuff heights • Combination of abutment with hex and castable sleeve
CLINICAL BENEFITS	<ul style="list-style-type: none"> • Screw retained abutments for crown or bridge restoration • A range of gingival heights • Castable sleeve for dentist/technician customization of the metal frame of the crown or bridge
MATERIAL	<p>Abutment: Titanium Alloy (Ti-6Al-4V ELI) Sleeve: Delrin</p>

ORDERING INFORMATION



Base	L (mm)	Screw	L (mm)
 NM-E1005	0.5	 NM-S1002	10.5
 NM-E1015	1.5	 NM-S1102	11.5
 NM-E1025	2.5	 NM-S1202	12.5

Hexed	Non-Hexed	L (mm)
 NM-C1015	 NM-C1007	12