


# CEMENT RETAINED RESTORATIONS |






## TEMPORARY ABUTMENTS



DESIGN FEATURES	<ul style="list-style-type: none"> <li>Sharp and deep retention grooves</li> <li>Available as hexed or non-hexed</li> </ul>
CLINICAL BENEFITS	<ul style="list-style-type: none"> <li>Acrylic retained abutment</li> <li>High acrylic retention.</li> <li>For single crown or bridge restoration</li> </ul>
MATERIAL	Titanium Alloy (Ti-6Al-4V ELI)

### ORDERING INFORMATION



Hexed	Non-Hexed	Short Non-Hexed	L (mm)	L1 (mm)	Ø D (mm)
		-	9.5	1.7	4.4
NM-A5001	NM-A5002				
	-		9.5	1.7	3.8
NM-A5101		NM-A5102			
	Included with all abutments above & available separately				
NM-S8324					