


OVERDENTURE RESTORATIONS |

FLAT ATTACHMENT



DESIGN FEATURES	<ul style="list-style-type: none"> • Flat attachment • Range of heights. • Components – metal house, wide range of nylon caps
CLINICAL BENEFITS	<ul style="list-style-type: none"> • Attachments for overdenture restoration • Range of gingival heights • Range of retention strength
MATERIAL	Titanium Alloy (Ti-6Al-4V ELI)

ORDERING INFORMATION



Flat Attachment	L (mm)	Flat Attachment	L (mm)
	0.5		3
NM-T1100		NM-T1103	
	1		4
NM-T1101		NM-T1104	
	2		5
NM-T1102		NM-T1105	

Set NM-T3099

Caps for Flat Attachment						
Metal	Extra Soft Nylon	Soft Nylon	Standard Nylon	High Retention Nylon	Laboratory	Sealing Disc
						
NM-T3010	NM-T3015	NM-T3016	NM-T3017	NM-T3018	NM-T3019	NM-T1824