

SCREW RETAINED RESTORATIONS |

MULTI-UNIT COMPONENTS

Healing Caps



Healing Caps
Ø 4.9mm H 4.8mm
NM-H7101



Healing Caps
Ø 4.9mm H 7.0mm
NM-H7102




Screw
NM-S7102

Transfers & Analog



Open Tray
Ø 4.9mm H 11mm
NM-T7111



Screw*
L 13mm
NM-S7111



Snap-On Cap for Ball Transfer
Ø 4.8mm H 11mm
NM-T4402



Closed Tray
Ø 4.9mm H 8mm
NM-T7102



Screw*
NM-S7102



Analog
Ø 4.9mm
NM-T7151

*Screw - Included with the transfer above & available separately


Abutments



Universal Abutment
Ø 4.9mm H 12mm
NM-T7121



Wide Universal Abutment
Ø 4.9mm H 12mm
NM-T7123



Castable Abutment
Ø 4.9mm H 12mm
NM-C7121



Screw
NM-S7102

Passive Fit Castable Abutment Set

The Passive Fit Castable Abutment Set consists of three parts aimed for the fabrications of accurate metal reinforced prostheses.

The Castable Sleeve Positioner is used for aligning the Castable Sleeve on the plaster model, ensuring passive fit of the fabricated metal cast when cemented to the Titanium Abutments.



Universal Abutment
Ø 4.9mm H 12mm
NM-T7121



Castable Sleeve
Ø 4.9mm H 10.5mm
NM-C7120



Castable Sleeve Positioner
Ø 4.9mm H 10mm
NM-T7122



Screw
NM-S7102

CAD-CAM COMPONENTS FOR MULTI UNIT

Titanium Bases



Ø 4.9mm H 4.5mm
NM-C7124



Ø 5.2mm H 4.5mm
NM-C7126



Screw
NM-S7102

Digital Analog for Multi-Unit



H 7mm
NMMDT7151

Scan Bodies



Ø 4.9mm H 7mm
NM-C9207



Ø 4.9mm H 10mm
NM-C9210



Ø 4.9mm H 13mm
NM-C9213



Screw
NM-S7102