

# SCREW RETAINED RESTORATIONS |

## Straight Multi-Unit



DESIGN FEATURES	<ul style="list-style-type: none"> <li>• Wide range of heights</li> </ul>
CLINICAL BENEFITS	<ul style="list-style-type: none"> <li>• Screw retained abutments for non-removable restorations</li> <li>• Wide range of gingival heights</li> </ul>
MATERIAL	Titanium Alloy (Ti-6Al-4V ELI)

### ORDERING INFORMATION

Straight Multi-Units					
Multi-Unit Platform Height (mm)					
1	2	3	4	5	
NM-A7101	NM-A7102	NM-A7103	NM-A7104	NM-A7105	



Abutment Carry-on handle  
NM-X7100

Included with all multi-unit bases

