

SCREW RETAINED & OVERDENTURE RESTORATIONS | VARI-CONNECT



DESIGN FEATURES	<ul style="list-style-type: none"> Two-piece angled multi-unit or angled attachment Wide range of angled multi-unit bases: 17°, 30°, 45°, 52° and 60° Wide range of multi-unit adapters heights. Ball-connect and flat-connect adapters Gold color anodization for soft tissue adaptation
CLINICAL BENEFITS	<ul style="list-style-type: none"> Angled screw retained abutments for non-removable restoration Angled attachments for overdenture restoration Wide range of gingival heights Great option for angled implants
MATERIAL	Titanium Alloy (Ti-6Al-4V ELI)

ORDERING INFORMATION

Angle Vari-Connect

Vari-Connect Platform Height (mm)

Angle°	2	3
--------	---	---

17°



NM-A7212



NM-A7213

30°



NM-A7233



NM-A7234

45°



NM-A7244

52°



NM-A7252

60°



NM-A7260

Vari-Connect Platform MULTI-UNIT Adapter

Multi-Unit Platform Adapter Height (mm)

	2	3	4	5	6
--	---	---	---	---	---



NM-T7282



NM-T7283



NM-T7284



NM-T7285



NM-T7286

Platform Adapter

